



INODAYA Hospitals - Kakinada

Documentation code:

INH/FMS.Doc.No:17

Policy on maintenance for piped medical gas, compressed air and vacuum installation

Prepared date:11/11/2025

Reference: FMS.6.f. NABH Standards – 6th Edition

Issue Date:11/11/2025

Issue no: 01

Review No: 0

Review date: 11/11/2026

PURPOSE:

The goal of the Hospital Standard is to ensure patient safety by continuous monitoring and improvement of safety aspects to deliver safe, effective and efficient care.

SCOPE:

The Hospital is committed to provide a safe Hospital environment for its External and Internal Customers.

RESPONSIBILITY:


BME Dept.

PLAN FOR COMPRESSED AIR:

- On a daily basis the Air compressor, Vacuum pump, Air filter, Electrical control panel and release drain line of air compressor is being checked.
- Preventive maintenance is being done every six months.
- Air sample calibration is being done TWICE in a year.

PLAN FOR VACUUM PUMP:

- On a daily basis vacuum pump, time V/S cycle and electrical panel is being checked.
- Mobil is replaced every three months.
- Preventive maintenance is done every six months.

Prepared by: 	Verified by: 	Approved by: 
Ms.Chandrakala	Mrs.G..Lakshmi Lavanya	Dr.G.Rammohan
Maintenance Incharge	Accreditation Coordinator	Managing Director